
Alternative Bowie-Dick-Test
PCD Class 2

EBI 16

STEAM

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According to EN 285
In compliance with EN ISO 11140-4

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Table of Contents

Safety Instructions	5
Unpacking/content of delivery	7
Overview	9
EBI16 checking the Vacuum test	11
Starting conditions.....	11
Place in the chamber.....	11
Evaluation	11
Archiving of the results	11
EBI16 checking the B&D test	13
Starting conditions.....	13
Place in the chamber.....	13
Evaluation	13
Archiving of the results	13
EBI16 checking the production cycle control	15
Starting conditions.....	15
Place in the chamber.....	15
Evaluation	15
Archiving of the results	15
EBI16 checking the special cycle control	17
Starting conditions.....	17
Place in the chamber.....	17
Evaluation	17
Archiving of the results	17
Required software	19
Logger and interface	19
Cool down times	19
Replacing the Battery	21
By using the software	21
Without software control	23
Replacement O-ring	25
What to do, if ...	27
Cleaning and Maintenance	29
Calibration service	29
Disposal	29
Technical Data	31
Approvals	31

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Safety Instructions

Do not use the data logger in areas in which there is a danger of explosions.

Danger to life!

Do not introduce the EBI16 into an environment that is hotter than 150°C.

The lithium battery may explode!

Do not subject the device to any microwave radiation.

The lithium battery may explode!

Do not grasp the EBI16 with bare hands immediately after use! Danger of burning!

You must only use the data logger with batteries supplied by **ebro**. These are the only batteries that have adequate heat resistance!

Data loggers emit wireless signals in the 2.4 GHz frequency range. These may cause significant interference to other technical devices! For this reason, only use data loggers in areas where these wireless emissions do not represent a risk!

Data loggers are currently only approved for use in Europe, USA, Canada and Japan.

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Unpacking/content of delivery

Check the device packaging and contents to ensure they are complete and undamaged. You must also check that the contents of the packaging match your order.

What is included in the delivery?

The delivery includes these parts:

- EBI16
- This operating manual

The delivery may also include various accessories.

- Programming interface (Standard) EBI IF150
- Programming interface EBI IF100/200
- Power supply AL 120
- USB cable
- Programming software Winlog.med
Winlog.validation
- Battery exchange set AL 103
- Rod antenna AL 111/112

If you have any reason for complaint, do not hesitate to contact us. You will find our contact details on the back of this manual

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Overview

The EBI16 data logger is a high precision device used to measure and record temperature and pressure.

It is powered by a lithium battery. The programming and analysis will be done with a PC

The EBI16 has no on/off switch. It is ready to be programmed until the battery runs out of power.

To program the EBI16, you will need a suitable interface (accessory) along with the appropriate read/write program on the PC (e. g. “*Winlog.med*”, also available as an accessory).

For more information about how to program the data logger and to extract the data, please refer to the user instructions for your read/write program.

The interface is connected to the PC via a universal serial bus (USB) port.

Always let the EBI16 cool down at room temperature.

The EBI16 has an internal cycle counter to indicate how many of the maximum 500 cycles have been run.

For processes with a high temperature phase > 30min, the electronics will be stressed more and so two cycles are counted. The operating time will be reduced accordingly.

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EBI16 checking the Vacuum test

Starting conditions

- The start temperature is not important for vacuum test.

Place in the chamber

- Place the EBI16 in the center of the Chamber.
- Place the EBI16 vertical in the chamber.
- Make sure that the EBI16 is not in the area of the steam entry.
- Start the sterilizer immediately after placement of the EBI16
- During the vacuum test, no other load is allowed in the chamber.

Evaluation

- Interface and software are required for the operation of the EBI16
- Automatic evaluation of the vacuum test can only be carried out with Winlog.med or Winlog.validation.

Archiving of the results

- The archiving of the results is ensured permanently by Winlog.med or Winlog.validation.

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
EBI16 checking the B&D test

Starting conditions

- EBI16 must be started at room temperature (+15°C ... +35°C).

Place in the chamber

- If necessary, preheat the sterilizer (sterilizer manual)
- Place the EBI16 in the center of the chamber.
- Place the EBI16 vertical in the chamber.
- Make sure that the EBI16 is not in the area of the steam entry.
- Start the sterilizer immediately after placement of the EBI16
- During the B&D test, no other load is allowed in the chamber.

Evaluation

- Interface and software are required for the operation of the EBI16
- Automatic evaluation of the B&D test can only be carried out with Winlog.med or Winlog.validation.

Archiving of the results

- The archiving of the results is ensured permanently by Winlog.med or Winlog.validation.

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EBI16 checking the special cycle control

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Place in the chamber

- If it is necessary, start before the actual sterilization program a vacuum test and/or a B&D test.
- Place the EBI16 in the Chamber. The EBI16 must be free without any obstruction. Do not place it inside a container or similar.
- Place the EBI16 vertical in the chamber.
- Make sure that the EBI16 is not in the area of the steam entry.
- Start the sterilizer immediately after placement of the EBI16.

Evaluation

- Interface and software are required for the operation of the EBI16
- Automatic evaluation of the special cycle control can be done only with Winlog.med or Winlog.validation.

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For processes with a high temperature phase > 30min, the electronics will be stressed extremely. Thereby, two cycles are counted. The operating time will be reduced accordingly

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
- Winlog.med V 3.5
- Winlog.validation V 3.5

Logger and interface

- To program or to start recording the EBI16 must be placed on the interface.
- To read out the data of the EBI16 it must be placed on the interface.
- Make sure that the case of the EBI16 is not hotter than +70°C. It could damage the housing of the interface otherwise

Cool down times

- To prevent damage to the housing of the interface, the EBI16 must cool down for approx. 7 minutes before it can be placed on the interface.
- To reach the start condition of +15°C ... +35°C -for programming- after a cycle with 134°C the EBI16 must cool down at room temperature for approx. 45 to 60 minutes.

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Replacing the Battery

An empty battery must be replaced.

Whilst changing the battery the operator must protect themselves against electrostatic discharges using suitable precautionary measures such as an earthed wristband!

Electrostatic discharges may destroy the device!

No data is lost when you change a battery or if the battery fails. However, if a power failure occurs, the built-in clock will stop and no longer show the correct time.

In this situation, extract the data from the data logger and then reprogram it. This synchronizes the date and time.

Recording continues immediately after you insert a new battery.

After changing the battery it must be activated by Software Winlog.med or Winlog.validation.

By using the software

Preferred option

For details in relation to the exchange of an empty battery for the EBI16 please consult the user manual of your evaluation software.

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



Without software control

To change a battery:


- Open the protection case of the EBI16. Data logger is now visible (see picture 1).
- Use the opening tool to open the housing of the data logger (see picture 2).
- Remove the battery retaining clip (picture 3). (Loosen the screw, lift the retaining clip and remove it carefully.)
- Pull the battery plug out of the circuit board (picture 4).
- Remove the empty battery.
- Place a new battery in the battery compartment in exactly the same position as the one you have just removed.
- Insert the retaining clip carefully in the two slots so that the red battery wire lies below it, so that it cannot be crushed or damaged.
- Plug the battery plug carefully into the socket on the circuit board (reverse battery protection). To the right of the battery compartment, you will see an LED. This LED will flash as soon as the circuit board is receiving power from the battery again.
- Tighten the retaining screw again.
- After the replacement of the O-ring reassemble the housing in reverse order.

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Replacement O-ring

- Remove the elastic silicone O-ring from the red plastic part of the logger housing (see picture 5).
- Cleaning the slot in which the O-ring was located (see picture 5).
- Take a new silicone O-ring, cover it with a thin layer of grease and place it in the slot.

The layer of grease should cover the O-ring completely (so that it glistens) but not leave lumps of residue on it.

A used silicone O-ring must always be disposed of. Under no circumstances should it be reused! This is because it is deformed when the two parts of the housing are screwed together. If it is used again it will not make a proper seal!

Take care not to overstretch the new elastic ring and not to damage when replacing it!

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What to do, if ...

- ... the data logger cannot be contacted from the interface:
 - Check that the interface has been connected correctly and that the evaluation software has recognized it.
 - Open the data logger (see also the pictures in the “Replacing the Battery” section).
 - Check to see whether condensation has penetrated the data logger housing.
 - Either wipe the condensation from them housing or allow it to air dry by leaving the data logger housing open.
 - Pull the battery plug from the socket and wait 30 seconds.
 - Then reconnect the battery and check that the LED to the right of the battery compartment is flashing. If the LED does not flash, the battery may be flat and needs to be replaced (see “Replacing the Battery”)
 - As a test, insert a fresh battery in the socket on the circuit board.
 - **If the LED still does not flash (despite the fresh battery), the problem lies with the data logger.**
 - Before you close the data logger housing again, replace the silicone O-ring.

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Cleaning and Maintenance

Use a damp cloth to clean the EBI16.

Do not use any solvents, such as acetone, because this damages the plastic.

Calibration service

The electronic B&D test is able for 500 cycles.

Please send the device for maintenance or calibration to our Service Department.

You will find our contact details on the back of this manual.

We shall return the fully calibrated device to you as soon as possible.

Disposal

If the device becomes no longer fit for purpose, it must be disposed of in a suitable, environmentally friendly manner.

Do not, under any circumstances, simply dispose of it in domestic garbage. You must return it to the manufacturer.

Dispose of the used battery in an environmentally friendly way at a suitable recycling point!

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Measure interval: 1 sec

Sensor:

Pressure piezo-resistive pressure sensor

Temperature PT 1000

Interface:

Wireless 2,4GHz / IEEE 802.15.4

Operating temperature 0 to +150°C

Storage -40 to +125°C

IP-rating IP 68 / NEMA 6P

Battery: Lithium 3.6 V, changeable

Battery lifespan:

up to 2 years dependent on application

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Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and

(2) This device must accept any interference received including interference that may cause undesired operation.

This device has been designed to operate with an internal antenna.

Other Antennas are strictly prohibited for use with this device.

Changes or modifications of the device could void the user's authority to operate the equipment.

Users manual or instruction manual may be downloaded at www.ebro.com over the internet.

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**Hersteller / Producer / Fabricant
& Service**

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WTW GmbH
WTW Wissenschaftlich-Technische Werkstätten GmbH